Temp Silicon Brush

Pennelli in silicone per modellazione composito
Silicon brushes for composite application
Silikoninstrumente für die Komposit-Modellierung
Pinceles de silicona para modelaciones en composite
Pinceau en silicone pour l'application du composite

MICERIUM SPA
Via G. Marconi, 83 - 16036 Avegno (Ge) Italy
Tel. 0185 7887 870

TPENN TPENS

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GB) Silicon brushes for composite application

- They do not stick on composite.
- They do not leave porosity and can be handled as normal brushes used by the lab.
- They have a resilience similar to uncured composite: this avoids air inclusion, and allows each layer of material to be compacted and smoothed.
- They can be easily cleaned: rinse excess soil from instrument. Disinfectant solution may be used in accordance with label instructions. Clean the instruments with a disposable tissue and disinfectant solution (as FD322 – DURR DENTAL), and leave them in an antibacterial disinfectant solution (as Peroxil or Vircon) for 10 minutes. Sold unsterilized.
- They can be used also with other materials: acrylics, monomers, sealers, lubricants, cianoacrylates, separators, gels, glazures, stains.

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INTRODUCTORY KIT TEMP BRUSHES black small rigid

5 assorted brushes Temp with silicon tip for the modellation of composite

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INTRODUCTORY KIT TEMP BRUSHES white small soft

5 assorted brushes Temp with silicon tip for the modellation of composite

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Silicon tip short spatula shape for modelling occlusal surfaces

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Silicon tip long spatula shape for blending layers and shaping surfaces

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Silicon tip round-flat shape for modelling cervical areas

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Silicon tip cone shape point ideal for occlusal shaping

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Silicon tip with three sides for stains and characterizations

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